



Pre-processing Testing Specification on Sugi® Material

| Tests/ items | Nominal values |
|---|--|
| Composition | |
| Cotton | Approx. 35 % |
| Regenerated cellulose | Approx. 65 % |
| Latex | None |
| Appearance/ Description | |
| Compressed white-cream absorbent sponge material | |
| Purity testing | |
| According to DAB or EP „Verbandswatte aus Baumwolle und Viskose“ in the current version. | |
| Surfactants | Foam must not completely cover the surface after 5 minutes |
| Water soluble substances | Max. 0.6% |
| Sulphide | Test solution must be coloured weaker than the reference solution |
| Drying loss | Max. 8% |
| Sulphate residue test (sulphate ash) | Max. 0.3% |
| Fluorescence | Light brown fluorescence and some yellow particles, however must not show strong blue, except for individual fibres. |
| Foreign fibres | none |
| Technical properties | |
| pH value of test solution | 6,5 – 7,5 |
| Water absorption | > 1200 % |
| Area weight range | Approx. 120 – 670g/m ² |
| Bioburden | |
| Bioburden limit for absorbent materials | < 200 cfu /5 g |
| Sterilization | |
| Ethylene oxide | Sterilization by ethylene oxide is recommended for Sugi Products. One-time re-sterilization with ETO does not impair the characteristics of the product. |
| Gamma rays | One-time re-sterilization with gamma rays will lead to a reduction of the absorptive capacity |
| Storage / Shelf life | |
| In sealed climate-secure foils at room temperature. Unsealed storage in a humid environment will lead to an increase of the dry layer thickness (absorption of humidity). | |
| Shelf life is 5 years. | |

ADDRESS

QUESTALPHA GmbH & Co. KG
Im Heerfeld 7
35713 Eschenburg
Deutschland/Germany

CONTACT

+49 (0) 2774 705 505
info@questalpha.com
www.questalpha.com

REGISTER/VAT

Amtsgericht Wetzlar
HRA 7726
VAT Number:
DE326299428

BANK DETAILS

Commerzbank AG
IBAN: DE42 5164 0043 0169 2540 00
BIC: COBADEFFXXX

GENERAL PARTNER

QUESTALPHA Verwaltungs-
GmbH, Eschenburg
Amtsgericht Wetzlar HRB 7958

CHAIRMAN

Jens T.W. Kuhn

MANAGING DIRECTOR

Dr. Sebastian Breuer

Biological Tests on Sugi® material

Cytotoxicity (Project 9111501, 1991; 97z021, 1997)

| | |
|----------------------------|--|
| Method: | ISO 10993-5 Biological evaluation of medical devices, Tests for <i>in vitro</i> Cytotoxicity. |
| Test Material: | SUGI® strips Dimension: 35 x 7.5 x 1.5 mm |
| Extract Production: | Extraction medium: phosphate buffered physiological saline solution, pH 7.4. Temperature and duration: 37 °C, 24 hours (1991) 37 °C, 72 hours (1997) |
| Cell Culture: | L 929 mouse fibroblasts |
| Controls: | Positive: Dilution series Negative: Extraction medium |
| Test Procedure: | 6 parallel cultures per tested dilution of the extracts (100, 60, 30, 10, 3, 1 and 0.3 %). Dyed with crystal violet. Calculation of cell growth inhibition from the extinction values at 580 nm |
| Results: | Test sample and negative control showed no signs of reactivity (Grade 0) with 30 % concentration of dilution |
| Conclusion: | SUGI® (non-sterilized and sterilized) does not exhibit pronounced cytotoxicity |

Sensitization (Report: 10-05-0923/00-92, 1992)

| | |
|----------------------------|--|
| Method: | ISO 10993-10 Biological evaluation of medical devices, Tests for irritation and skin sensitization. |
| Test Material: | SUGI® strips Dimension: 35 x 7.5 x 1.5 mm |
| Extract Production: | Extraction medium: phosphate buffered 0.9% saline solution, pH 7.4. Temperature and duration: 37 °C, 72 hours |
| Test Population: | 20 guinea pigs |
| Control Population: | 10 guinea pigs |
| Test Procedure: | Induction exposure to undiluted extract of test material. After Extract production with surface to volume ratio 3 cm ² /ml. Release 14 days afterwards, with this extract |
| Results: | After 24, 48 and 72 hours no allergic reaction, no toxic effects. Classification numbers for erythema and oedema = 0. |
| Conclusion: | SUGI® is non-sensitizing |

Irritation – Skin (Report: 10-03-0922/01-92, 1992)

| | |
|-------------------------|--|
| Method: | ISO 10993-10 Biological evaluation of medical devices, Tests for irritation and skin sensitization. |
| Test Material: | Undiluted extract made of SUGI® strips |
| Test Population: | 3 rabbits |
| Test procedure: | duration of contact: 4 hours |
| Results: | After 30 minutes, 1, 24, 48 and 72 hours., no allergic reaction, no toxic effects. Irritation index for animals = 0. |
| Conclusion: | SUGI® is non-irritating to the skin |

Irritation – Ocular (Report: 10-03-0937/00-92, 1992)

| | |
|-------------------------|--|
| Method: | ISO 10993-10 Biological evaluation of medical devices, Tests for irritation and skin sensitization. |
| Test Material: | Undiluted extract made of SUGI® strips |
| Test Population: | 3 rabbits |
| Test Procedure: | duration of contact: 4 hours |
| Results: | After 1, 24, 48 and 72 hours, no allergic reactions, no toxic effects or eye lesions. Classification number = 0. |
| Conclusion: | SUGI® is non-irritating to the ocular surface |

Update biocompatibility

Cytotoxicity (Project 14Z057, 2014)

| | |
|----------------------------|---|
| Method: | ISO 10993-5 Biological evaluation of medical devices, Tests for <i>in vitro</i> Cytotoxicity. |
| Test Material: | SUGI® macro swabs (34 mm, non-compressed), SUGI® sponge strips (35 x 7.5 x 1.5 mm, compressed) |
| Extract Production: | Sterile test items were transferred into the eluent (0.2 per milliliter cell culture medium containing 10 % fetal calf serum) in consideration of absorption capacities: SUGI® macro swabs 14.8 ml/g, SUGI® sponge strips 15 ml/g. Temperature and duration: 37 °C, 24 hours. |
| Cell Culture: | L 929 mouse fibroblasts |
| Controls: | Positive: Dilution series of Dimethylsulfoxide Negative: Extraction medium |
| Test Procedure: | 96 parallel cultures per tested dilution of the extracts (100, 30, 10, 3 %). After 4 hours incubation filled with dilution series and incubates for 72 hours at 37 °C again. Dyed with crystal violet. Calculation of cell growth inhibition from the extinction values at 570 nm |
| Results: | Test samples and negative controls showed no signs of reactivity (Grade 0) |
| Conclusion: | SUGI® does not exhibit pronounced cytotoxicity |

GC/MS Fingerprint (Project 14Y102, 2014)

| | |
|----------------------------------|--|
| Method: | ISO 10993-18 Biological evaluation of medical devices, Chemical characterization of materials. |
| Test Material: | SUGI® macro swabs (34 mm, non-compressed), SUGI® sponge strips (35 x 7.5 x 1.5 mm, compressed) |
| Sample Extraction: | Medium: water, isopropyl alcohol, n-hexane. Temperature and duration: 37 °C, 72 hours. In closed glass vials |
| Controls: | Positive: solution of n-tetradecane dissolved in n-hexane Negative: Pure extraction |
| Test Procedure: | 1 µl of each extract was injected into the GC (duplicate analysis per vial), separated on the capillary column and detected by MS (mass selective detector). |
| Results & Conclusion: | No semi-volatile organic compounds were detected above the analysis limit |