
Sugi®

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
Stand der Information 28.02.2024

Sugi®

Sterile absorbent material for use in the field of surgery (single use)


Sugi® Rhino
(sterile absorbent swabs with retrieval cord according to Prof. Dr. H. Stammberger)

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

Sugi® Endopad (sterile)

Ø 34 mm
according to Prof. Dr. H. Stammberger
(6 sachets à 1)
REF 31785

Sugi® OtoSwab
with retrieval cord (sterile)

Ø 3 x 1 mm, expanded 3 x 5 mm

60 swabs (10 sachets à 6)
REF 31651

Rhino Swabs

Rhino Long

Endopad

OtoSwab

Composition

Cellulose
Regenerated cellulose
Cord: polyester

Properties

- Sugi absorbent materials are manufactured from high-quality cellulose material and are toxicologically harmless.
- The swabs/pads expand immediately on contact with liquid and become soft and pliable.
- Sugi absorbent materials exhibit high (> 50 N) tear strength.
- Sugi absorbent materials will not adhere to wounds if they are removed within 24 hours.
- Sugi absorbent materials leave a clean wound surface after removal.
- The absorptive capacity is approximately 17 times its own weight.

a) Rhino absorbent swabs

Indication

Conditions of the nose requiring interventions

Range of application

- Preoperative and intraoperative nasal swabs for use in nasal surgery especially during:
 - Septoplasty
 - Septorhinoplasty
 - Conchotomy
 - Paranasal sinus surgery
- Swabs for use in endoscopic ENT surgery
- Swabs for use in skull base surgery involving nasal surgery procedures
- Carrier medium for liquid antibiotic or vasoconstrictor medication

b) Endopads

Indication

- Intraoperative cleaning of optical instruments (e.g. endoscopes) and simultaneous application of an anti-fogging agent.

c) OtoSwab

Indication

Conditions of the ear requiring interventions

Range of application

- Intra-operative swab for ear surgery
 - Postauricular approach to middle ear surgery
 - Endaural approach to middle ear surgery
 - Meatoplasty
 - Tympanoplasty
 - Tympanotomy
 - Cortical mastoidectomy
 - Middle ear surgical procedures
 - Cochlear Implant
- Carrier medium for vasoconstrictive/antibiotic solutions

Contraindications and Incompatibility

None known

Application

a) Rhino absorbent swabs */OtoSwab

The swabs are soaked in a vasoconstrictor and/or anesthetic solution for preparing the operation site. They are then placed on the relevant mucosal areas for 5 – 10 minutes.

They are easily removed using the attached retrieval cord.

The swabs can also be used intraoperatively to absorb blood and secretion.

Manufacturer:
QUESTALPHA GmbH & Co. KG
Im Heerfeld 7
35713 Eschenburg, Germany
www.questalpha.com

QUESTALPHA
For Medical and Life Science

Distributed in the U.S. by:
Kettenbach LP
62-64 Enter Lane
Islandia, NY 11749, USA

Made in Germany
37248/0924

CE 0123

b) Endopads

The dry pad is prepared by soaking it in an anti-fogging agent. After removing the protective foil from the adhesive strip, the pad is positioned close to the operation site.

Caution!

- Sugi absorbent products must be disposed of in the same way as dressing material and must not be reused, as there is a risk of infection and the function and safety of the product is no longer guaranteed.
- Prior to use, visually inspect the sterile barrier system for damage and do not use it if damaged or already opened.
- The Sugi absorbent materials shown are radiolucent.
- **Sugi absorbent materials are to be used by medical personnel qualified to perform or assist in surgery.**
- * In rare instances the physicochemical conditions associated with nasal surgery, both with and without nasal packing, may present a risk of toxic shock syndrome (TSS). Warning signs of TSS include: sudden high fever, vomiting, diarrhea, dizziness, fainting (or near fainting when standing up), and/or a rash that looks like a sunburn. Should one of these symptoms occur, administer the appropriate medical treatment.
- The container of the Sugi Endopad must not be glued to skin.

Intended patient group

Patients requiring the absorption of body fluids or topical administration of liquids

Legend

Keep away from sunlight	
Use-by date	
Sterilized using ethylene oxide	STERILE EO
Do not re-use	
Keep dry	
Batch Code	LOT
Dual system (Germany only)	
Catalogue number	REF
Consult instructions for use	
Do not use if package is damaged	
Temperature limit	
Medical device	MD
Date of manufacture	

Limitation of liability

Except where prohibited by law, QUESTALPHA GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the legal basis of the claim, including warranty, contract, negligence, or strict liability.

The information provided for QUESTALPHA products is based on comprehensive research and experience in application technology.

Results are furnished to the best of our knowledge, subject to technical changes within the framework of product development.

However, users must comply with and consider all recommendations and information in connection with any use.

In the event of serious incidents in connection with the product, please contact the manufacturer and the responsible authority immediately.

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