
# Sugi®

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Distributed in the U.S. by:  
Kettenbach LP  
62-64 Enter Lane  
Islandia, NY 11749, USA

Made in Germany  
37248/1423

QUESTALPHA  
For Medical and Life Science

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
Stand der Information 2. Februar 2023

# Sugi®

Sterile absorbent material for use in the field of surgery (single use)


**Sugi® Rhino**  
(sterile absorbent swabs with retrieval cord according to Prof. Dr. H. Stammberger)

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

### OtoSwab

#### Composition

Cellulose  
Regenerated cellulose  
Cord: polyester

#### Properties

- Sugi absorbent materials are manufactured from high-quality cellulose material and are toxicologically harmless.
- The swabs/pads expand immediately on contact with liquid and become soft and pliable.
- Sugi absorbent materials exhibit high (> 50 N) tear strength.
- Sugi absorbent materials do not adhere to the wound surface and are therefore easily removed.
- Sugi absorbent materials leave a clean wound surface after removal.
- The absorptive capacity is approximately 17 times its own weight.

#### Range of Applications

##### a) Rhino absorbent swabs

- Preoperative and intraoperative nasal swabs for use in nasal surgery especially during:
  - Septoplasty
  - Septorhinoplasty
  - Conchotomy
  - Paranasal sinus surgery
- Swabs for use in endoscopic ENT surgery
- Swabs for use in skull base surgery involving nasal surgery procedures
- Carrier medium for liquid antibiotic or vasoconstrictor medication

##### b) Endopads

- Intraoperative cleaning of optical instruments (e.g. endoscopes) and simultaneous application of an anti-fogging agent.

##### c) OtoSwab

- Intra-operative swab for ear surgery
  - Postauricular approach to middle ear surgery
  - Endaural approach to middle ear surgery
  - Meatoplasty
  - Tympanoplasty
  - Tympanotomy
  - Cortical mastoidectomy
  - Middle ear surgical procedures
  - Cochlear Implant
- Carrier medium for vasoconstrictive/antibiotic solutions

#### Contraindications and Incompatibility

None known

#### Application

##### a) Rhino absorbent swabs \*/OtoSwab

The swabs are soaked in a vasoconstrictor and/or anesthetic solution for preparing the operation site. They are then placed on the relevant mucosal areas for 5 – 10 minutes.

They are easily removed using the attached retrieval cord.

The swabs can also be used intraoperatively to absorb blood and secretion.

##### b) Endopads

The dry pad is prepared by soaking it in an anti-fogging agent. After removing the protective foil from the adhesive strip, the pad is positioned close to the operation site.

Manufacturer:  
**QUESTALPHA GmbH & Co. KG**  
Im Heerfeld 7  
35713 Eschenburg, Germany  
www.questalpha.com

**QUESTALPHA**  
For Medical and Life Science

Distributed in the U.S. by:  
**Kettenbach LP**  
62-64 Enter Lane  
Islandia, NY 11749, USA

Made in Germany  
37248/1423

#### Caution!

- Sugi absorbent materials are to be disposed of like bandages and must not be reused, as there is a risk of infection and the function and safety of the product are no longer guaranteed.
- Do not use sterile products if the sterile packaging is damaged or has already been opened.
- The Sugi absorbent materials shown are radiolucent.
- Sugi absorbent materials are to be used by medical personnel qualified to perform or assist in surgery.**
- \* In rare instances the physiochemical conditions associated with nasal surgery, both with and without nasal packing, may present a risk of toxic shock syndrome (TSS). Warning signs of TSS include: sudden high fever, vomiting, diarrhea, dizziness, fainting (or near fainting when standing up), and/or a rash that looks like a sunburn. Should one of these symptoms occur, administer the appropriate medical treatment.
- The container of the Sugi Endopad must not be glued to skin.

#### Intended patient group

Patients requiring the absorption of body fluids or topical administration of liquids

#### Legend

|                                  |            |
|----------------------------------|------------|
| Keep away from sunlight          |            |
| Use-by date                      |            |
| Sterilized using ethylene oxide  | STERILE EO |
| Do not re-use                    |            |
| Keep dry                         |            |
| Batch Code                       | LOT        |
| Dual system (Germany only)       |            |
| Catalogue number                 | REF        |
| Consult instructions for use     |            |
| Do not use if package is damaged |            |
| Temperature limit                |            |
| Medical device                   | MD         |
| Date of manufacture              |            |

#### Limitation of liability

Except where prohibited by law, QUESTALPHA GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the legal basis of the claim, including warranty, contract, negligence, or strict liability.

The information provided for QUESTALPHA products is based on comprehensive research and experience in application technology.

Results are furnished to the best of our knowledge, subject to technical changes within the framework of product development.

However, users must comply with and consider all recommendations and information in connection with any use.

In the event of serious incidents in connection with the product, please contact the manufacturer and the responsible authority immediately.

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